


HUMANE SLAUGHTER OF SHEEP



Electrode position
(front view)



Electrode position
(side view)



Stick position
1 (recommended)
2 (alternative)

MINIMUM STUNNING CURRENT: 1 amp adult & 0.6 amp lambs

ELECTRODE APPLICATION TIME: NOT LESS than 3 seconds

STUN-TO-STICK TIME: LESS than 15 seconds

Phase	Duration	Visible signs	Action
TONIC	10 to 20 sec	<ul style="list-style-type: none"> > animal collapses > body becomes rigid > no rhythmic breathing > ears raised > hind legs flexed into body 	<ul style="list-style-type: none"> > stick during this phase if possible
CLONIC	1 to 45 sec	<ul style="list-style-type: none"> > involuntary kicking or paddling > relaxation 	<ul style="list-style-type: none"> > stick immediately
RECOVERY	30 to 60 sec	<ul style="list-style-type: none"> > resumption of normal rhythmic breathing > responds to painful stimuli > becomes visually aware > attempts to stand 	<ul style="list-style-type: none"> > stun with captive-bolt and stick immediately

For further information contact:

Humane Slaughter Association

The Old School, Brewhouse Hill, Wheathampstead, Herts AL4 8AN, UK

Tel: +44 (0)1582 831919 Fax: +44 (0)1582 831414 E-mail: info@hsa.org.uk

www.hsa.org.uk. Charity Registered in England No. 209563 © HSA 2008